

GRAM

SERIES / SERIE / SÉRIE

RZ


EN

ES

FR

CE

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INDEX

PRECAUTIONS	5
BEFORE USING THE BALANCE	5
PRODUCT INTRODUCTION	6
• Features & Specifications	6
• Front panel	7
• Keyboard & functions of keys	7
POWER SUPPLY	9
SETTING PARAMETERS	9
CALIBRATION	12
OPERATION	13
• Weighing	13
• Tare	13
• Checking function	14
• Counting pieces	15
• Accumulation	15
CONNECT THE BALANCE TO A PRINTER	16
EXPORTING DATA TO PC IN FORM EXCEL	17
U.KEY CONNECTOR	19
GARANTIE	20

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PRECAUTIONS

- * Place the scale on a flat and stable surface.
- * Verify that the input voltage and the plug type matches the local AC power supply.
- * Warm up the scale for 15 minutes before using it for the first time.
- * Keep the scale away from EMI noise, strong wind and vibration, which might cause incorrect reading.
- * Avoid sudden temperature changes (suitable operating temperature is between 0°C~ 40°C.)
- * Disconnect the power supply when cleaning the scale.
- * Do not immerse the scale in water or other liquids.
- * Service should be performed by authorized personnel only.

BEFORE USING THE BALANCE

Unpacking and Checking

Open the package and check the instrument for transport damage. Immediately inform your dealer if you have complaints or if parts are missing. The package should contain:

- Scale body
- Weighing pans
- User manual
- Adaptor

Installing Components

1) Before using the scale, remove the shipping protection screw (rotate counterclockwise), which located underneath the scale and cork the plug buckle.

Note: the RZ-30 model doesn't include the shipping protection screw.

2) Cover the weighing pans on the scale body properly.


Leveling the Scale

To compensate for small irregularities or inclinations at the location, the scale can be leveled. The scale is equipped with a level indicator at the front panel .Adjust the leveling feet until the air bubble in the indicator is centered as shown.

Note: The scale should be leveled each time once its location is changed.

PRODUCT INTRODUCTION

1. Leveler
2. Keypad
3. Adjusting feet
4. Plastic weighing pan
5. Stainless Steel weighing pan
6. On/Off switch
7. Transport protection screws
8. Power socket


SPECIFICATIONS


FEATURES

- Resolution up to 1/15000
- Tare, Preset Tare, auto zero tracking
- Single point and linear calibration available
- HI/OK/LO quantity checking function
- Weighing platter easy to disassemble with dust proof PVC cover
- Low battery and charging status indication
- Accumulation of weight and quantity
- RS 232 port

Model		RZ-6	RZ-15	RZ-30
Capacity(kg)		6	15	30
Weighing units	Kg/g/lb			
Display	LCD backlit displays (6/5/6 Digits)			
Pan size	294 x 228 x 13.5mm			
Scale dimensions	330 x 289 x 104mm			
Power supply	DC 9v/400mA or rechargeable battery (6V/4A)			

FRONT PANEL

-  Low battery indicator
-  Tare indicator
-  Zero indicator. The zero range is $\pm 2\%$ maximum capacity
-  Battery charge indicator
-  High limit indicator (HI)
-  Ok limit indicator.
-  Low limit indicator (LOW).
-  Predetermined Tare symbol “▼”
-  Net weight symbol “▼” it light on when a tare is done or a predetermined tare value has been entered.
-  Stable symbol .
-  “M+” symbol “▼” when weighing accumulation is activated.
-  Unit weight symbol, it light on when an unknow unit weight has been entered

KEYBOARD

Short press steps through activated weighing units, release on desired one.

Displays gross and net weight by turns.

1. Tares the weight of the container or accepts the keypad tare entries.
2. Preset Tare and cancel the Tare
3. Select the later parameter in the same level
4. Select the later record when checking the accumulation record
5. Shift keys rightwards

1. Zeros the display (within 2% of max. capacity)
2. Exits from certain operation without save

1. Long press to enter function setting
2. Deletes accumulation records
3. Delete the present accumulation record for the sake of wrong accumulation operation under the accumulation mode

1. Adds the indicated weight into accumulation memory
2. During editing, save and return to the higher option

1. Select the former parameter in the same level
2. Select the former record when checking the accumulation record
3. Shift keys leftwards.

1. During setting value, add 1 to the current value
2. Set the upper limit of the check weighing
3. Print

1. During setting value, deduct 1 to the current value
2. Set the lower limit of the check weighing
3. Hold function

POWER SUPPLY

Please verify the local power source before plugging into the power outlet, and use the individual power socket and original adaptor.

Alternative power supply

- 1) (9V/400mA) adaptor
- 2) (6V/4A) internal rechargeable battery

Power Consumption

About 300mW (without backlight)

About 380 mW (with backlight)

Low battery warning

When “” appears in the upper left corner of the weight window, the battery power requires recharging. The charge lamp turns green from red when the recharging is completed (which takes about 8 hours). Disconnect the scale from power supply when it is fully charged. A full charged battery can operate the scale for approximately 65 hrs with backlight, 80hrs without backlight.

Note: Battery is to be replaced only by an authorized service dealer .Risk of explosion can occur if replaced with the wrong type or connected improperly

SETTING MODE

Function Setting

- 1) Press  key while powering on, or  key in the weighing mode to enter parameter setting. The window displays “check”
- 2) Press key  or  to shift between the functions
- 3) Press key  to enter the selected parameter.
- 4) Press key  or  to circle the items of present parameter.
- 5) Press key  to save after changing and return to parameters choosing.
- 6) Press key  to return without saving.
- 7) Press key  again to return to weighing mode.

PARAMETERS

P00 *CHEFF* Internal code checking

P01 *LIGHT* Backlight modes

on = Backlight

off = No backlight

Auto = Auto-on after getting a stable weigh data or any key is pressed. But auto off after N seconds (N=2s, 5s, 10s, 20s, ever) with no action.

(ever=the backlight lasts forever with more than 9d placed on the pan) when it is stable.

P02 *Unit* Current using units setting

Unit: Press the  key to select the default unit when powering the balance: "pcs, lb, oz, g. kg. final" Final= keep the final being used unit when power off

Use: Press the  key to select the weighing unit

On = Enable the unit

Off = Disable the unit

P03 *PERO* Zero range setting

d0 (one division) , d1 (2 divisions), d2 (3 divisions), d3 (4 divisions), d4 (5 divisions) and d5 (6 divisions)

P04 *Hold* Hold function

Hold 0 : No hold function

Hold 1 : Peak hold

Hold 2: Hold after stable

Hold 3: Hold after stable. Release after moving away the article. The hold value is base on the current value and its range could be set in sub menu. Accumulation hold function is available, that is you could add article after hold the first value.

Hold 4: Press the key 

Sub-menu for Hold 3 mode: INF (default: infinity) / 10 / 20 / 50 / 100 / 200 / 500 / 1000 / 2000/5000/10000/20000/50000

H=current hold value, R=hold value range, d=division, W=actual weight

P05 *Stable* Stable weight checking setting

On: it starts checking weight when the weight and quantity are within the checking range and the stable symbol appears.

Off: it starts checking weight once the weight and quantity are within the checking range

P06 *Check* Check weighing function

On: Check weighing under the condition that the weight is within the limits and the stable indication appears.

Off: Check weighing under the condition that the weight is within the limits.

P07 *Beep* Check weighing buzzer beep

Hi: There will be a warning sound when the weight of articles exceeds the upper limit and the weight is equal or more than 20d.

Lo: There will be a warning sound when the weight of articles exceeds the lower limit and the weight is equal or more than 20d

Out: There will be a warning sound when the weight of articles is between the upper and lower limit and the weight is equal or more than 20d

No beep: No beep

P08 *PC* External devices

PC = Computer

U.KEY: U.KEY connector works with PC to send the weighing data to computer (WinXP/Win7) in form of Excel and so on.

P09 *Baud* Baud rate setting: 9600、4800、2400

P10 *Print* Printing modes setting

contin: Printing continually

stable: Printing automatically when get a stable weight (more than 20d)

key: Printing by pressing keys

P11 *Print* Printing format setting: prt01~prt03

P12 *Filter* Filtering setting: 1, 2, 3, 4 Set the filtering level in which the stable indication

turns on .The higher the setting, the slower stabilization time.

P13 Precondition for zeroing or taring setting:

stable : Taring or zeroing can be not action without stable indication after pressing the key  or 

always: Taring or zeroing can be action without stable indication after pressing the key  or 

auto: Press key  or  , but it works only when get the stable indication.

P14 Weight memory parameter

On: display the last weight when powering on again

Off: do not display the last weight when powering on again

P15 Parameter initialization

Press key  then key  to start initialization. It is finished when displaying OK.

CALIBRATION

Note: Before calibration, please set the capacity first. The unit used in calibration is the one that has been set before. During the calibration procedure, long press  to return to normal weighing mode without saving.

Here we take 3kg/0.2g as an example

1. Press and hold  while powering on, not release it till the window displays “**CAL**”
2. With no load on the weighing pan, press  to start zero point calibration.
3. Wait till the window displays the first calibration value.

Note: The first calibration value is default. If the value has been changed, the default value is 1/3 of full load.

Note: If you need to change the value, do as the following: Press key  to enter the value setting. Press key  or  to move leftwards or rightwards. Press key  or  to modify the value (or via numeric keys). Press key  to save.

4. Put the corresponding weight on the weighing pan, and then press  to complete the first point calibration.

Note: After the first point calibration, the window can display the weight value. If no need for the other point calibration, press the  to save and exit from calibration mode.

5. Add another weight 500g to the current weight. The window will show the total weights value, press the  key to complete

Repeat this step to complete the calibrations of on3, on4 and on5.

6. Press  to save. After the window displays "PASS", It will return to normal weighing mode.

Note: the window displays  automatically when on5 is finished. It will be saved and return to the normal weighing.

OPERATION

Weighing

Put item to be weighed on the scale. The Weight window shows the weight.

Tare

When weighing a sample that must be held in a container, tare and store the container weight into memory.

1) Under the weighing mode, put the container on the scale.

Press the key  complete tare action.

The  symbol appears and the window displays as shown on the right

2) Put the item into the container, then the window displays the net weight.

3) Clear the pan. The window displays the taring weight, that is the weight of container.

4) To cancel tare action, clear the pan and press  or  keys

Preset Tare

1) Long press the  key for three seconds, it enters to preset tare setting.

Preset-tare value setting: Press key  or  to shift and key  or  to modify the value (or via numeric keys).

Press key  to save and return to weighing mode.

2) Put the items into the container and the weight of container will be tared automatically from the total weight.

NOTE: Press the key  to display the gross and net weight in turns.

3) To cancel the preset-tare action, clear the pan and press  or  keys

Checking Function

Upper weight limit

1) Long press key  under the weighing mode.

The number leftmost is blinking.

2) Press key  or  to shift and key  or  to modify the value

3) Press key  to confirm and save.

Lower limit weight

1) Long press key  under the weighing mode.

The number leftmost is blinking.

2) Press key  or  to shift and key  or  to modify the value

3) Press key  start or close weighing checking

4) Press key  confirm and save.

Put the sample on the weighing pan.

HI indication appears, when the item on the weighing pan is greater than the upper limit

OK indication appears, when the item on the weighing pan is between upper and lower limits.

LO indication appears, when the item on the weighing pan is smaller than lower limit

Note: the item on the weighing pan should be more than or equal to 20e.

When changing the Hi-Lo value, the scale will activate the weighing checking function automatically. If the Lo value is higher than Hi value, then the Hi value will become the same value as Lo data.

Basic counting pieces

- 1) Under the weighing mode, press the  key to select the unit “PCS”
- 2) Long press the  key, the ex-factory default sample quantity (10pcs) is displayed.
- 3) Use the  and  keys to choose the sampling amount. Available options are: 10, 20, 50, 100, 200, 500, 1000 (pieces)
- 4) Put the corresponding samples on the weighing pan and then press the  key. The message  is displayed momentarily before the display reverts to the sample quantity.
- 5) Remove the samples and put the load on, the scale calculates the amount of the load.
- 6) To go back to the normal weighing mode, remove the load and press the  key to select the proper unit weight.

Accumulation

- 1) Under the weighing mode, put the item on the weighing pan. Press the  key, the display will show momentarily the message “Add01”.
- 2) Remove the item and the display goes back to zero before the next accumulation can register. (The maximum is 99 pieces, display the lastest 10 accumulation events in details)

Accumulation display

- 1) Press key  to displays the total accumulation data.
- 2) Press key  or  check the total accumulation or each accumulation in detail.

The number of accumulated weight is showed at the left bottom.

Note: XX=total accumulation events; XX=the single accumulation

Accumulation clear

To clear the total accumulation data or the first 10 accumulation events, press key  to display the accumulation data needed and then press the key  to clear data chosen.

When clear total accumulation data, the accumulation signal “▼” will disappear and back

to normal weighing mode. If need to exit without clearing data and return to normal weighing mode, press the  key.

INITIALIZATION OF THE SCALE TO PRINTER (Optional)

1. Press the  key while powering on to enter parameter setting. Use the  and  to shift to P09 parameter.

2. Press  key to enter the external device setting. Press  or  keys to choose the printer model.

3. After choosing one printer model, press the  key and the display will show “UNSUP” or “INIT?”

- “UNSUP” means the printer can’t be initialized by scale. Press the  key to return.
- “INIT” means the printer is should be initialized. Press the  key to initialize the printer. When the initialization is finished, the display will show the message “ok” and then displays the printer model. Press the  key to return.

Input commands

Connect the scale to computer. Run serial port debugging program on the computer. Input the capital number “Z”, “T”, “R”, “C”, “P” in the sending area and the indicator can conduct the corresponding actions:

Z = Zero

T = Tare


C = Cancel Tare

R/P = Reading / Print

SERIAL INTERFACE

If external interface is needed, please select proper RS-232 board and only when this board is adopted, the functions can be enabled.

Note: RTC Is unavailable.

Troubleshooting and Error message

Error Message	Problem	Solution
ERR0	Exceed the zero range	The item on the pan should be within 2% of full load.
ERR2	Initial zero point exceeds 30% of full load	1.Remove the obstacle article under the pan 2.Replace the load cell or contact the maintenance department.
ERR3	Exceed the A/D resolution range	1. Replace A/D 2..Replace the load cell or contact the maintenance department.
ERR4	EEPROM failure	Re-weld EEPROM or contact the maintenance department.
ERR5	Overload condition	Remove the overload item
ERR6	Exceeds the display range	-----
ERR7	Accumulated number of weighments exceeds the display range	Delete the exceeding weights
ERR8	Weight limit value is higher than the full load value	Reset the weight limit value.
ERR9	Exceed tare or pre-tare range	$0 < \text{Tare value} \leq \text{full load}$
ERR10	Wrong calibration weights	Place the right test weights and the calibration value should be below full load.

Exporting data to PC in the form of EXCEL

Introduction:

Connect the scale with PC and set the parameter of external device as "EXCEL" on the scale, then you could export the weighing data to PC in the form of EXCEL. With this function, you could record/accumulate/average/data statistical analysis the testing data, which we could call it as **scale-computer data management function**.

Note: pls enable "Use Serial Keys" function in the computer.


Hardware connection and settings

1. Use transmitting serial wire or USB wire to connect scale and pc.
Note: pls install usb driver first, if you use usb wire.
2. Parameter settings in scale:
"PERI" = "EXCEL" (external device)
"BAUD" = "2400"/"4800"/"9600" (baud rate)
"PRT.M" = "KEY"/"STABLE" (printing model)
"PRT.F" = "PRT.F01" (printing format)


Enable the function of “Use Serial Keys” in the computer

Set Windows XP as a example:

1. Press “Start” ->“Run”, and enter “access.cpl” ->“OK”.


2. In the dialog box “Accessibility Options”, enter General option, choose “Use Serial Keys” and press “Settings”.
In the dialog box “Settings for SerialKeys”, set the corresponding Serial port and Baud rate, which should be same as Baud rate in scale.

3. Test if Serial Keys works well.
Open a Text Document, and press the Print button on the scale. The Serial Keys works well, if pc exports the weighing data to Text.

Export weighing data to Excel

1. Open Excel.
2. Press [Print], then Excel will show Weighing data.
3. Use “Format Cells” to beautify Excel:

Troubleshooting

1. Enable “Fast User Switching” function will disable “Use Serial Keys” in the computer.
Note: Start->"Control Panel"->"User Account"-> “Fast User Switching” function.
2. The data export requires only “GND+RX+TX” three lines, while some wire include nine lines (including usb 9-25pin adapter), which may cause abnormal function.
Note: cut other useless line to check if these three lines work well.

U.KEY Connector

Introduction

Connect the scale with PC and set the parameter of external device as “U.KEY” on the scale, then you could export the weighing data to PC in the form of EXCEL. With this function, you could record/accumulate/average/data statistical analysis the testing data, which we could call it as **scale-computer data management function**.

Most importantly, you don't need any setup in computer, but only use U.KEY connector to link computer system directly, such as EXCEL, Express System, ERP system and so on. This connector is compatible with windows2000/xp/Vista/win7/win8/Linux/Android system, and the weighing data will show on the position where the mouse cursor in.

With one USB port, you could connect scale to computer easily.

Precaution

1. Close the input method program before connecting.
2. Ensure the decimal point digit is same as the parameter setting of scale, if use EXCEL file

in computer.

For Example:

Scale shows: Total Weight: 1.000KG 、 Unit Weight: 0.001g、 Total: 1000PCS

EXCEL Setting: the decimal point digit should be 3 for total weight and unit weight column, and the decimal point digit should be 0 for total column.

Setting

1. Prepare one computer with windows2000/xp/Vista/win7/win8/Linux system or Android pad.

2. Scale setting:

"PERI" = "U.KEY" (External devices)

"BAUD" = "9600" (Baud)

"PRT.M" = "KEY"/"STABLE" (Print Mode)

"RTC" = "ON"/"OFF" (Time)

3. Open the receiving port such as EXCEL or ECR system, and the weighing data will show on the position where the mouse cursor in

WARRANTY

This scale is guaranteed for one year from the delivery date. The guarantee covers any fabrication defect of the material.

During this period GRAM PRECISION, covers the manpower and the spare parts necessary for the reparation of the scale.

This guarantee does not cover the failures caused by an inappropriate use or overcharge.

The guarantee does not cover the freight cost (transport) necessary to repair the scale.

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- © Pour obtenir une précision optimale, placer la balance sur une surface plane, ferme et sans vibration.